



Product Name: Recombinant Human OLR1 Protein
Catalog #: PHH3116



Summary

Name	OLR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P78380
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQVSDLLTQE QANLTHQKKK LEGQISARQQ AEEASQSEEN ELKEMIETLA RKLNEKSKEQ MELHHQNLNL

image



Background

Alternative Names	Oxidized low-density lipoprotein receptor 1 (Ox-LDL receptor 1) (C-type lectin domain family 8 member A) (Lectin-like oxidized LDL receptor 1) (LOX-1) (Lectin-like oxLDL receptor 1) (hLOX-1) (Lectin-type oxidized LDL receptor 1) [Cleaved into: Oxidized low-density lipoprotein receptor 1, soluble form]
--------------------------	--

Note

Product Name: Recombinant Human OLR1 Protein
Catalog #: PHH3116



For Research Use Only , Not for Diagnostic Use