

Product Name: Recombinant Human ROR1 Protein
Catalog #: PHH3127



Summary

Name	ROR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01973
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLLALLAAL LLAARGAAQ ETELSVSAEL VPTSSWNISS ELNKDSYLT DEPMNITTS LGQTAELHCK VSGNPPPTIR WFKNDAPVVQ EPRRLSFRST IYGSRLRIRN LDTTDTGYFQ CVATNGKEVV SSTGVLFVKF GPPPTASPGY SDEYEEDGFC QPYRGIACAR FIGNRTVYME SLHMQGEIEN QITAAFTMIG TSSHLSDKCS QFAIPSLCHY AFPYCDTSS VPKPRDLCDR ECEILENVLC QTEYIFARSN PMILMRLKLP NCEDLPQPES PEAANCIRIG IPMADPINKN HKCYNSTGVD YRGTVSVTKS GRQCQPWNSQ YPHHTHTFTAL RFPENGGHS YCRNPGNQKE APWCFTLDEN FKSDLCDIPA CDSKDSKEKN

image



Background

Alternative Names Inactive tyrosine-protein kinase transmembrane receptor ROR1

Product Name: Recombinant Human ROR1 Protein
Catalog #: PHH3127



(Neurotrophic tyrosine kinase, receptor-related 1)

Note

For Research Use Only , Not for Diagnostic Use