



Product Name: Recombinant Human DSG1 Protein
Catalog #: PHH3130



Summary

Name	DSG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02413
Host	HEK293
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TYVVTGNMGS NDKVGDFVAT DLDTGRPSTT VRYVMGNNPA DLLAVDSRTG KLTLKNKGTK EQYNMLGGKY QGTILSIDDN LQRTCTGTIN INIQSFGNDD RTNTEPNTKI TTNTGRQUEST SSTNYDTSTT STDSSQVYSS EPGNGAKDLL

image



Background

Alternative Names	Desmoglein-1 (Cadherin family member 4) (Desmosomal glycoprotein 1) (DG1) (DGI) (Pemphigus foliaceus antigen)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use