



Product Name: Recombinant Human POU4F2 Protein
Catalog #: PHH3147



Summary

Name	PO4F2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q12837
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDADPRDLEA FAERFKQRRRI KLGVTQADVG SALANLKIPG VGSLSQSTIC RFESLTLSHN NMIALKPILQ AWLEEEAEKSH REKLTKEPELF NGAEKKRKRT SIAAPEKRSL EAYFAIQPRP SSEKIAAIAE KLDLKKNVVR VWFCNQRQKQ

image



Background

Alternative Names POU domain, class 4, transcription factor 2 (Brain-specific homeobox/POU domain protein 3B) (Brain-3B) (Brn-3B)

Note

For Research Use Only , Not for Diagnostic Use