


Product Name: Recombinant Human PTPRJ Protein
Catalog #: PHH3155



Summary

Name	PTPRJ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12913
Host	HEK293
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPIPDPSVAT VATGENGITQ ISSTAESFHK QNGTGTPQVE TNTSEDGESS GANDSLRTPE QGSNGTDGAS QKTPSSTGPS PVFDIKAVSI SPTNVILTWK SNDTAASEYK YVVKHKMENE KTITVVHQPW CNITGLRPAT SYVFSITPGI GNETWGDPRV IKVITEPIPV SDLRVALTGV RKAALSWSNG NGTASCRVLL ESIGSHEELT QDSRLQVNIS GLKPGVQYNI NPYLLQSNKT KGDPLGTEGG LDASNTERSR AGSPTAPVHD ESLVGPVDPS SGQQSRDTEV LLVGLEPGTR YNATVYSQAA NGTEGQPQAI EFRTNAIQVF DVTAVNISAT SLTLIWKVSD NESSSNITYK IHVAGETDSS NLNVSEPRAV IPGLRSSTFY NITVCPVLGD IEGTPGFLQV HTPPVPSDF

image



Background

Product Name: Recombinant Human PTPRJ Protein
Catalog #: PHH3155



Alternative Names

Receptor-type tyrosine-protein phosphatase eta (Protein-tyrosine phosphatase eta) (R-PTP-eta) (EC 3.1.3.48) (Density-enhanced phosphatase 1) (DEP-1) (HPTP eta) (Protein-tyrosine phosphatase receptor type J) (R-PTP-J) (CD antigen CD148)

Note

For Research Use Only , Not for Diagnostic Use