



Product Name: Recombinant Human GRIN2B protein
Catalog #: PHH3165



Summary

Name	NMDE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q13224
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDAHEKDDFH HLSVVPVEL VAMNETDPKS IITRICDLMS DRKIQGVVFA DDTDQEIAIQ ILDFISAAQL TPILGIHGGG SMIMADKDES

image



Background

Alternative Names	Glutamate receptor ionotropic, NMDA 2B (GluN2B) (Glutamate [NMDA] receptor subunit epsilon-2) (N-methyl D-aspartate receptor subtype 2B) (NMDAR2B) (NR2B) (N-methyl-D-aspartate receptor subunit 3) (NR3) (hNR3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use