



Product Name: Recombinant Human MSLN Protein
Catalog #: PHH3170



Summary

Name	MSLN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q13421
Host	HEK293
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGQEAAPLDG VLANPPNISS LSPRQLLGFP CAEVSGLSTE RVRELAVALA QKNVKLSTEQ LRCLAHRLSE PPEDLDALPL DLLLFLNPDA FSGPQACTRF FSRITKANVD LLPRGAPERQ RLLPAALACW GVRGSLLSEA DVRALGGLAC DLPGRFVAES AEVLLPRLVS CPGPLDQDQQ EAARAALQGG GPPYGPPSTW SVSTMDALRG LLPVLGQPPII RSIPQGIVAA WRQRSSRDPS

image



Background

Alternative Names	Mesothelin (CAK1 antigen) (Pre-pro-megakaryocyte-potentiating factor) [Cleaved into: Megakaryocyte-potentiating factor (MPF); Mesothelin, cleaved form]
--------------------------	--

Product Name: Recombinant Human MSLN Protein
Catalog #: PHH3170



Note

For Research Use Only , Not for Diagnostic Use