


Product Name: Recombinant Human SLC44A2 Protein
Catalog #: PHH3220



Summary

Name	CTL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8IWA5
Host	HEK293
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKVIYPTDSR GEFCGQKGTK NENKPYLFYF NIVKCASPLV LLEFQCPTPQ ICVEKCPDRY LTYLNARSSR DFEYYKQFCV PGFKNNKGVA EVLQDGDGCPA VLIPSKPLAR RCFPAIHAYK GVLMVGNETT YEDGHGSRKN ITDLVEGAKK ANGVLEARQL AMRIFEDYTV

image



Background

Alternative Names Choline transporter-like protein 2 (Solute carrier family 44 member 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC44A2 Protein
Catalog #: PHH3220

