

## Summary

<b>Name</b>	BMPER_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8N8U9
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SFLTGSVAKC ENEGEVLQIP FITDNPCIMC VCLNKEVTCK REKCPVLSRD CALAIKQRGA CCEQCKGCTY EGNTYNSSFK WQSPAEPCLV RQCQEGVVTE SGVRCVVHCK NPLEHLGMCC PTCPGCVFEG VQYQEGEEFQ PEGSKCTKCS CTGGRTQCVR EVCPII LSCPQ HLSHIPPGQC CPKCLGQRKV FDLPGSCLF RSDVYDNGSS FLYDNCTACT CRDSTVVCKR KCSHPGGCDQ

## image



## Background

<b>Alternative Names</b>	BMP-binding endothelial regulator protein (Bone morphogenetic protein-binding endothelial cell precursor-derived regulator) (Protein crossveinless-2) (hCV2)
--------------------------	--

**Product Name: Recombinant Human BMPER Protein**  
**Catalog #: PHH3230**



---

**Note**

For Research Use Only , Not for Diagnostic Use