



Product Name: Recombinant Human PRG4 Protein
Catalog #: PHH3242



Summary

Name	PRG4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q92954
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDGLTTLRNG TLVAFRGHYF WMLSPFSPPS PARRITEVWG IPSPIDTVFT RCNCEGKTFE FKDSQYWRFT NDIKDAGYPK PIFKGFGLT GQIVAALSTA KYKNWPESVY FFKRGGSIQQ YIYKQEPVQK

image



Background

Alternative Names Proteoglycan 4 (Lubricin) (Megakaryocyte-stimulating factor) (Superficial zone proteoglycan) [Cleaved into: Proteoglycan 4 C-terminal part]

Note

For Research Use Only , Not for Diagnostic Use