



Product Name: Recombinant Human IL17F Protein
Catalog #: PHH3248



Summary

Name	IL17F_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PD4
Host	HEK293
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTVKT LHGPA MVKYL LLSIL GLAFLSEAAA RKIPKVGHTF FQKPESCPPV PGGSMKLDIG IINENQRVSM SRNIESRSTS PWNYYTVTWDP NRYPSEVVQA QCRNLGCINA QGKEDISMNS VPIQQETLVV RRRKHQGCSVS FQLEKVLVTV GCTCVTPVIH HVQ

image



Background

Alternative Names Interleukin-17F (IL-17F) (Cytokine ML-1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IL17F Protein
Catalog #: PHH3248

