


Product Name: Recombinant Human CD274 Protein
Catalog #: PHH3272



Summary

Name	PD1L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His & C-HA
Accession #	Q9NZQ7
Host	HEK293
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTVPKDLYVV EYGSNMTIEC KFPVEKQLDL AALIVYWEME DKNIIQFVHG EEDLKVQHSS YRQRARLLKD QLSLGNAALQ ITDVKLQDAG VYRCMISYGG ADYKRITVKV NAPYNKINQR ILVVDPTSE HELTCQAEY PKAEVIWTSS DHQVLSGKTT TTNSKREEKL FNVSTLRIN TTTNEIFYCT FRRLDPEENH TAEVIPELP

image



SDS



WB

Background

Alternative Names	Programmed cell death 1 ligand 1 (PD-L1) (PDCD1 ligand 1) (Programmed death ligand 1) (hPD-L1) (B7 homolog 1) (B7-H1) (CD antigen CD274)
--------------------------	--

Note

Product Name: Recombinant Human CD274 Protein
Catalog #: PHH3272



For Research Use Only , Not for Diagnostic Use