



Product Name: Recombinant Human KLK13 Protein
Catalog #: PHH3284



Summary

Name	KLK13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKR3
Host	HEK293
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGVSQESSKV LNTNGTSGFL PGGYTCTPHS QPWQAALLVQ GRLLCGGVLV HPKWWLTAAH CLKEGLKVYL GKHALGRVEA GEQVREVVHS IPHPEYRRSP THLNHDHDIM LLELQSPVQL TGYIQTLPIS HNNRLTPGTT CRVSGWGTTT SPQVNYPKTL QCANIQLRSD EECRQVYPGK ITDNMLCAGT KEGGKDSCEG DSGGPLVCNR TLYGIVSWGD FPCGQPDRPG VYTRVSRYVL WIRETIRKYE TQQQKWLGKP Q

image



Background

Alternative Names	Kallikrein-13 (EC 3.4.21.-) (Kallikrein-like protein 4) (KLK-L4)
--------------------------	--

Product Name: Recombinant Human KLK13 Protein
Catalog #: PHH3284



Note

For Research Use Only , Not for Diagnostic Use