


Product Name: Recombinant Human TRPM3 Protein
Catalog #: PHH3328



Summary

Name	TRPM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9HCF6
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPYMSQAQEI HLQEKEAEEP EKPTKEKEEE DMELTAMLGR NNGESSRKKD EEEVQSKHRL IPLGRKIYEF

image



Background

Alternative Names	Transient receptor potential cation channel subfamily M member 3 (Long transient receptor potential channel 3) (LTrpC-3) (LTrpC3) (Melastatin-2) (MLSN2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use