


**Product Name: Recombinant Human CD4 Protein**  
**Catalog #: PHH3355**



## Summary

<b>Name</b>	CD4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P01730
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	42 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LLVLQLALLP AATQGKKVVL GKKGDTVELT CTASQKKSQ FHWKNSNQIK ILGNQGSFLT KGPSKLNDR A DSRRSLWDQG NFPLIKNLK IEDSDTYICE VEDQKEEVQL LVFGLTANS D THLLQGQSLT LTLESPPGSS PSVQCRSPRG KNIQGGKTLS VSQLELQDSG TWTCTVLQNQ KKVEFKIDIV VLAFAQASSI VYKKEGEQVE FSFPLAFTVE KLTGSGELWW QAERASSSKS WITFDLKNKE VSVKRVTQDP KLQMGKKLPL HLTLPQALPQ YAGSGNLTLA LEAKTGKLHQ EVNLVVMRAT QLQKNLTCEV WGPTSPKMLL SLKLENKEAK VSKREKAVVW LNPEAGMWQC LLSDSGQVLL ESNIKVLPTW

## image



## Background

**Alternative Names** T-cell surface glycoprotein CD4 (T-cell surface antigen T4/Leu-3) (CD antigen)

**Product Name: Recombinant Human CD4 Protein**  
**Catalog #: PHH3355**



---

CD4)

**Note**

For Research Use Only , Not for Diagnostic Use