



**Product Name: Recombinant Human CD44 Protein**  
**Catalog #: PHH3376**



## Summary

<b>Name</b>	CD44_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P16070
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MGRGHQAGRR MDMDSSHSIT LQPTANPNTG LVEDLDRTGP LSMTTQQSNS QSFSTSHEGL EEDKDHPTTS TLTSSNRNDV TGRRDPNHS EGSTTLLEGY TSHYPHTKES

## image



## Background

<b>Alternative Names</b>	CD44 antigen (CDw44) (Epican) (Extracellular matrix receptor III) (ECMR-III) (GP90 lymphocyte homing/adhesion receptor) (HUTCH-I) (Heparan sulfate proteoglycan) (Hermes antigen) (Hyaluronate receptor) (Phagocytic glycoprotein 1) (PGP-1) (Phagocytic glycoprotein I) (PGP-I) (CD antigen CD44)
--------------------------	--

## Note

**Product Name: Recombinant Human CD44 Protein**  
**Catalog #: PHH3376**



---

For Research Use Only , Not for Diagnostic Use