


Product Name: Recombinant Human TNFRSF8 Protein
Catalog #: PHH3380



Summary

Name	TNR8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P28908
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGTIPQAKPT PVSPATSSAS TMPVRRGGTRL AQEAASKLTR APDSPSSVGR PSSDPGLSPT QPCPEGSGDC

image



Background

Alternative Names Tumor necrosis factor receptor superfamily member 8 (CD30L receptor) (Ki-1 antigen) (Lymphocyte activation antigen CD30) (CD antigen CD30)

Note

For Research Use Only , Not for Diagnostic Use