



Product Name: Recombinant Human MUC16 Protein
Catalog #: PHH3381



Summary

Name	MUC16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8WXI7
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQTSTQLPTT TSAHPGQVPR SAATTLDVIP HTAKTPDATF QRQGQTALTT EARATSDSWN EKEKSTPSAP WITEMMNSVS

image



Background

Alternative Names Mucin-16 (MUC-16) (Ovarian cancer-related tumor marker CA125) (CA-125) (Ovarian carcinoma antigen CA125)

Note

For Research Use Only , Not for Diagnostic Use