

Product Name: Recombinant Human CD99 Protein
Catalog #: PHH3382



Summary

Name	CD99_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14209
Host	HEK293
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALPDNENKKP TAIPKKPSAG DDFDLGDAVV DGENDDRPP NPPKMPNPN PNHPSSSGSF SDADLADGVS

image



Background

Alternative Names CD99 antigen (12E7) (E2 antigen) (Protein MIC2) (T-cell surface glycoprotein E2) (CD antigen CD99)

Note

For Research Use Only , Not for Diagnostic Use