


Product Name: Recombinant Human DCTN1 Protein
Catalog #: PHH3404



Summary

Name	DCTN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q14203
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EASARPLRVG SRVEVIGKGH RGTVAIVGAT LFATGKWVGV ILDEAKGKND GTVQGRKYFT CDEGHGIFVR QSQIQVFEDG ADTTSPETPD

image



SDS

Background

Alternative Names	Dynactin subunit 1 (150 kDa dynein-associated polypeptide) (DAP-150) (DP-150) (p135) (p150-glued)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use