

Product Name: Recombinant Human CDON Protein
Catalog #: PHH3439




Summary

Name	CDON_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4KMG0
Host	HEK293
Species	Human
Molecular Mass	102 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FTSEPLSAVQ KLGGPVVLHC SAQPVTTTRIS WLHNGKTLDG NLEHVKIHQG TLTILSLNSS LLGYYQCLAN NSIGAIVSGP ATVSVAVLGD FGSSTKHVIT AEEKSAGFIG CRVPESNPKA EVRYKIRGKW LEHSTENYLI LPSGNLQILN VSLEDKGSYK CAAYNPVTHQ LKVEPIGRKL LVSRPSSDDV HILHPTHSQA LAVLSRSPVT LECVVSGVPA PQVYWLKDGQ DIAPGSNWRR LYSHLATDSV DPADSGNYSC MAGNKSQDVK YVTYMNVLVLE HASISKGLQD QIVSLGATVH FTCDVHGNPA PNCTWFHNAQ PIHPSARHLT AGNGLKISGV TVEDVGMYYQC VADNGIGFMH STGRLEIEND GGFKPVIITA PVSQKADGD FVTLSCNASG LPVPVIRWYD SHGLITSHPS QVLRKSRKS QLSRPEGLNL EPVYFVLSQA GASSLHIQAV TQEHAGKYIC EAANEHGTQ AEASLMVVPF ETNTKAETVT LPDAAQNDDR SKRDGSETGL LSSFPVKVHP SAVESAPEKN ASGISVPDAP IILSPPQHTT PDTYNLVWRA GKDGGLPINA YFVKYRKLLD GVGMLGSWHT VRVPGSENEL HLALEPSSL YEVLMMVARSA AGEQGPAMLT FRTSKEKTAS SKNTQASSPP VGIPKYPVVS EAANNFNGVV LTDSSRHSGV PEAPDRPTIS TASETSVYVT WIPRANGGSP ITAFKVEYKR MRTSNWLVAE EDIPPSKLSV EVRSLPGST YKFRVIAINH YGESFRSSAS RPYQVVGFPN RFSSRPITGP HIAYTEAVSD TQIMLKWTYI PSSNNNTPIQ GFYIYRPTD SDNDSYKRD VVEGSKQWHM IGHLPQETSY DIKMQCFNEG GESEFSNMI CETKVCRVPG ASEYPVKDLS TPPNSLGSGG NVGPATSPAR

image

Product Name: Recombinant Human CDON Protein
Catalog #: PHH3439



SDS



WB

Background

Alternative Names

Cell adhesion molecule-related/down-regulated by oncogenes

Note

For Research Use Only , Not for Diagnostic Use