


Summary

Name	ANGP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15123
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEYDDSVQRL QVLENIMENN TQWLMKLENY IQDNMKKEMV EIQQNAVQNNQ TAVMIEIGTN LLNQTAEQTR KLTDVEAQVL NQTRLELQL LEHSLSTNKL EKQILDQTSE INKLQDKNSF LEKKVLAMED KHIIQLQSIK EEKDQLQVLV SKQNSIIEEL EKKIVTATVN NSVLQKQHD LMETVNNLLT MMSTSNSAKD PTVAKEEQIS FRDCAEVFKS GHITNGIYTL TFPNSTEEIK AYCDMEAGGG GWTIIQRRED GSVDFQRTWK EYKVGFGNPS GEYWLGNEFV SQLTNQQRVY LKIHLKDWEG NEAYSLEYHF YLSSEELNYR IHLKGLTGTA GKISSIONPG NDFSTKDGDN DKCICKCSQM LTGGWWFDAC GPSNLNGMYYPQRQNTNKFN GIKWYYWKGS GYSLKATTMM

image



Background

Product Name: Recombinant Human ANGPT2 Protein
Catalog #: PHH3481



Alternative Names

Angiopoietin-2 (ANG-2)

Note

For Research Use Only , Not for Diagnostic Use