



Product Name: Recombinant Human APOE3 protein
Catalog #: PHH3513



Summary

Name	APOE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P02649
Host	HEK293
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVEQAVETEP EPELRQQTEW QSGQRWELAL GRFWDYLRWV QTLSEQVQEE LLSSQVTQEL RALMDETMKE LKAYKSELEE QLTPVAEETR ARLSKELQAA QARLGADMED VCGRLVQYRG EVQAMLGQST EELRVRLASH LRKLRKRLLR DADDLQKRLA VYQAGAREGA ERGLSAIRER LGPLVEQGRV RAATVGSLAG QPLQERAQAW GERLRARMEE MGSRTDRDL EVKEQVAEVR AKLEEQAQOI RLQAEAFQAR LKSWFEPLVE DMQRQWAGLV EKVQAAVGTS AAPVPSDNH

image



Background

Alternative Names	Apolipoprotein E (Apo-E)
--------------------------	--------------------------

Product Name: Recombinant Human APOE3 protein
Catalog #: PHH3513



Note

For Research Use Only , Not for Diagnostic Use