


**Product Name: Recombinant Human ATP2B4 Protein**  
**Catalog #: PHH3543**



## Summary

<b>Name</b>	AT2B4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P23634
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	AIAKCKQAGI TVRMVTGDNI NTARAIATKC GILTPGDDFL CLEGKEFNRL IRNEKGEVEQ EKLDKIWPKL RVLARSSPTD KHTLVKGIID STVGEHRQVV AVTGDGTNDG PALKKADVGF

## image



## Background

<b>Alternative Names</b>	Plasma membrane calcium-transporting ATPase 4 (PMCA4) (EC 7.2.2.10) (Matrix-remodeling-associated protein 1) (Plasma membrane calcium ATPase isoform 4) (Plasma membrane calcium pump isoform 4)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ATP2B4 Protein**  
**Catalog #: PHH3543**

---

