



Product Name: Recombinant Human NRG1 Protein
Catalog #: PHH3593



Summary

Name	NRG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-hFc & C-His
Accession #	Q02297
Host	HEK293
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKGKKKERGS GKKPESAAGS QSPALPPRLK EMKSQESAAG SKLVLCETS SEYSSLRFKW FKNGNELNRK NKPQNIQK KPGKSELRLN KASLADSGEY MCKVISKLGN DSASANITIV ESNEITGMP ASTEGAYVSS ESPIRISVST EGANTSSSTS TSTTGTSHLV KCAEKEKTCF VNGGECFMVK DLSNPSRYLC KCQPGFTGAR CTENVPMKVQ NQEKAEELYQ

image



Background

Alternative Names	Pro-neuregulin-1, membrane-bound isoform (Pro-NRG1) [Cleaved into: Neuregulin-1 (Acetylcholine receptor-inducing activity) (ARIA) (Breast cancer cell differentiation factor p45) (Glial growth factor) (Heregulin) (HRG) (Neu differentiation factor) (Sensory and motor neuron-derived factor)]
--------------------------	---

Product Name: Recombinant Human NRG1 Protein
Catalog #: PHH3593



Note

For Research Use Only , Not for Diagnostic Use