


Product Name: Recombinant Human MUCL3 protein
Catalog #: PHH3646



Summary

Name	MUCL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q3MIW9
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHIVLHSGQR PPELPKSTEI HEQKRHCNTT RHSKPTDKPT GNSKTIDHKS STDNHEAPPT SEENSSNQGK DPMIRNQRSV DPADSTTTHK ESAGKKHITP

image



Background

Alternative Names Mucin-like protein 3 (Diffuse panbronchiolitis critical region protein 1)

Note

For Research Use Only , Not for Diagnostic Use