


Product Name: Recombinant Human VSTM1 Protein
Catalog #: PHH3662



Summary

Name	VSTM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UX27
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKNEKPPKPS LHAWPSSVVE AESNVTLCQ AHSQNVTFVL RKNVDSGYKQ EQSSAENEAE FPFTDLKPKD AGRYFCAYKT TASHEWSESS EHLQLVVTDK HDELEAPSMK

image



SDS



WB

Background

Alternative Names	V-set and transmembrane domain-containing protein 1 (Signal inhibitory receptor on leukocytes-1) (SIRL-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use