


**Product Name: Recombinant Human MILR1 Protein**  
**Catalog #: PHH3670**



## Summary

<b>Name</b>	MILR1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-hFc & C-His
<b>Accession #</b>	Q7Z6M3
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	51 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KAVLDCEAMK TNEFPSPCLD SKTKVVMKGQ NVSMFCSHKN KSLQITYSLF RRKTHLGTQD GKGEPAIFNL SITEAHESGP YKCKAQVTSC SKYSRDFSFT IVDPVTSPVL NIMVIQTETD RHITLHCLSV NGSLPINYTF FENHVAISPA ISKYDREPAE FNLTKKNPGE EEEYRCEAKN RLPNYATYSH PVTMPSTGGD

## image



## Background

**Alternative Names** Allergin-1 (Allergy inhibitory receptor 1) (Mast cell antigen 32) (MCA-32) (Mast cell immunoglobulin-like receptor 1)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human MILR1 Protein**  
**Catalog #: PHH3670**

---

