


Product Name: Recombinant Human CLSPN Protein
Catalog #: PHH3736



Summary

Name	CLSPN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9HAW4
Host	HEK293
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRLVSNDNEF DSDEDEHSDS GNDLALEDHE DDDEEELLKR SEKLKMQMRL RKYLEDEAEV SGSDVGSSEDE YDGEEIDEYE EDVIDEVLPS DEELQSQIKK IHMKTMLDDD KRQLRLYQER YLADGDLHSD GPGRMRKFRW KNIDDASQMD LFHRDSDDDQ TEEQLDESEA RWRKERIERE QWLRDMAQQG KITAEIEEEI

image



Background

Alternative Names Claspin (hClaspin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CLSPN Protein
Catalog #: PHH3736

