



Product Name: Recombinant Human EDAR Protein
Catalog #: PHH3770



Summary

| | |
|------------------------------|---|
| Name | EDAR_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9UNE0 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 17 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | TGLCQECPPC GPGEEPYLSC GYGTKDEDYG CVPCPAEKFS KGGYQICRRH KDCEGFFRAT VLTPGDMEND AECGPCLPGY YMLENRPRNI YGMVCYSCLL APPNTKECVG ATSGASANFP GTSGSSTLSP FQHAHKELSG |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Tumor necrosis factor receptor superfamily member EDAR (Anhidrotic ectodysplasin receptor 1) (Downless homolog) (EDA-A1 receptor) (Ectodermal dysplasia receptor) (Ectodysplasin-A receptor) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EDAR Protein
Catalog #: PHH3770

