


Product Name: Recombinant Human CSPG5 Protein
Catalog #: PHH3804



Summary

Name	CSPG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O95196
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEGETPSWSL LDLYDDFTPF DESDFYPTTS FYDDLDEEEE EEEDDKDAVG GGDLEDENEL LVPTGKPLGL PGTGQPTRSW HAVPPQHTLG

image



Background

Alternative Names Chondroitin sulfate proteoglycan 5 (Acidic leucine-rich EGF-like domain-containing brain protein) (Neuroglycan C)

Note

For Research Use Only , Not for Diagnostic Use