


Product Name: Recombinant Human F3 Protein
Catalog #: PHH3814



Summary

Name	TF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P13726
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTDEIVKDVK QTYLARVFSY PAGNVESTGS AGEPLYENSP EFTPYLETNL GQPTIQSFEQ VGTKVNVTV

image



Background

Alternative Names	Tissue factor (TF) (Coagulation factor III) (Thromboplastin) (CD antigen CD142)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use