



Product Name: Recombinant Human IFNLR1 Protein
Catalog #: PHH3828



Summary

Name	INLR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8IU57
Host	HEK293
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTWLPGLGNP QDVTFVAYQ SSPTRRRWRE VEECAGTKEL LCSMMCLKKQ DLYNKFKGRV RTVSPSSKSP WVESEYLDYL FEVEPAPPVL VLTQTEEILS ANATYQLPPC MPPLDLKYEY AFWKEGAGNK TLFVTPHGQ PVQITLQPAA SEHHCLSART IYTFVSPKYS KFSKPTCFL

image



Background

Alternative Names	Interferon lambda receptor 1 (IFN-lambda receptor 1) (IFN-lambda-R1) (Cytokine receptor class-II member 12) (Cytokine receptor family 2 member 12) (CRF2-12) (Interleukin-28 receptor subunit alpha) (IL-28 receptor subunit alpha) (IL-28R-alpha) (IL-28RA) (Likely interleukin or cytokine receptor 2) (LICR2)
--------------------------	--

Product Name: Recombinant Human IFNLR1 Protein
Catalog #: PHH3828



Note

For Research Use Only , Not for Diagnostic Use