


Summary

Name	LRP12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y561
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGVSTACGET PEQIRAPSGI ITSPGWPSEY PAKINCSWFI RANPGEIITI SFQDFDIQGS RRCNLDWLT I ETYKNIESYR ACGSTIPPPY ISSQDHIWIR FHSDDNISRK GFRLAYFSGK SEEPNCACDQ FRCGNGKCIP EAWKCNNMDE CGDSSDEEIC AKEANPPTAA AFQPCAYNQF QCLSRFTKVY TCLPESLKCD GNIDCLDLGD EIDCDVPTCG QWLKYFYGTF NSPNYPDFYP PGSNCTWLID TGDHRKVILR FTDFKLDGTG YGDYVKIYDG LEENPHKLLR VLTAFDSHAP LTVVSSSGQI RVHFCADKVN AARGFNATYQ VDGFLPWEI PCGGNWGCYT EQQRCDGYWH CPNGRDETNC TMCQKEEFPC SRNGVCYPRS DRCNYQNHCP NGSDEKNCFF CQPGNFHCKN NRCVFESWVC DSQDDCGDGS DEENCPVIVP

image



Background

Product Name: Recombinant Human LRP12 Protein
Catalog #: PHH3835



Alternative Names

Low-density lipoprotein receptor-related protein 12 (LDLR-related protein 12) (LRP-12) (Suppressor of tumorigenicity 7 protein)

Note

For Research Use Only , Not for Diagnostic Use