


Product Name: Recombinant Human IFNGR1 Protein
Catalog #: PHH3856



Summary

Name	INGR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15260
Host	HEK293
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPTNVTIESY NMNPIVYWEY QIMPQVPVFT VEVKNYGVKN SEWIDACINI SHHYCNISDH VGDPSNSLWV RVKARVGQKE SAYAKSEefa VCRDGKIGPP KLDIRKEEKQ IMIDIFHPSV FVNGDEQEVD YDPETTCYIR VYNVYVRMNG SEIQYKILTQ KEDDCDEIQC QLAIPVSSLN SQYCVSAEGV LHVWGVVTEK SKEVCITIFN

image



Background

Alternative Names	Interferon gamma receptor 1 (IFN-gamma receptor 1) (IFN-gamma-R1) (CDw119) (Interferon gamma receptor alpha-chain) (IFN-gamma-R-alpha) (CD antigen CD119)
--------------------------	---

Product Name: Recombinant Human IFNGR1 Protein
Catalog #: PHH3856



Note

For Research Use Only , Not for Diagnostic Use