


Product Name: Recombinant Human NOTCH3 Protein
Catalog #: PHH3865



Summary

Name	NOTC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	tag free
Accession #	Q9UM47
Host	HEK293
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSEPRCPRA ACQAKRGDQR CDRECNSPGC GWDGGDCSLS VGDPWRQCEA LQCWRLFNNS RCDPACSSPA CLYDNFDCHA GGRERTCNPV YEKYCADHFA DGRCDQGCNT EECGWDGLDC ASEVPALLAR GVLVLTVLLP PEELLRSSAD FLQRLSAILR TSLRFRLDAH GQAMVFPYHR PSPGSEPRAR RELAPEVIGS VVMLEIDNRL CLQSPENDHC FPDAQSAADY LGALSAVERL DFPYPLRDVR GEPLEPPEPS

image



Background

Alternative Names	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]
--------------------------	---

Product Name: Recombinant Human NOTCH3 Protein
Catalog #: PHH3865



Note

For Research Use Only , Not for Diagnostic Use