



Product Name: Recombinant Human ATP1B4 Protein
Catalog #: PHH3866



Summary

Name	AT1B4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9UN42
Host	HEK293
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLNGFLQGYN DSLQEEMNVD CPPGQYFIQD GNEDEDKKAC QFKRSFLKNC SGLEDPTFGY STGQPCILLK MNRIVGFRPE LGDPVKVSCK VQRGDENDIR SISYYPESAS FDLRYPPYYG

image



Background

Alternative Names	Protein ATP1B4 (X,K-ATPase subunit beta-m) (X/potassium-transporting ATPase subunit beta-m)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use