


Product Name: Recombinant Human CD96 Protein
Catalog #: PHH3871



Summary

Name	TACT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P40200
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQPSNSSMTT RGFNYPWTSS GTDTKKSISR IPSETYSSSP SGAGSTLHDN VFTSTARAFS EVPTTANGST KTNHVHITGI

image



Background

Alternative Names	T-cell surface protein tactile (Cell surface antigen CD96) (T cell-activated increased late expression protein) (CD antigen CD96)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use