



Product Name: Recombinant Human KIAA0319L Protein
Catalog #: PHH3882



Summary

Name	K319L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZA0
Host	HEK293
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADAGPDKELT LPVDSTTLDG SKSSDDQKII SYLWEKTQGP DGVQLENANS SVATVTGLQV GTYVFTLTVK DERNLQSQSS VNVIVKEEIN KPPIAKITGN VVITLPTSTA ELDGSKSSDD KGIVSYLWTR DEGSPAAGEV LNHSDDHPIL FLSNLVEGTY TFHLKVTDK GESDTRTTV EVKPDPRKNN LVEILDINV

image



Background

Alternative Names	Dyslexia-associated protein KIAA0319-like protein (Adeno-associated virus receptor) (AAVR)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KIAA0319L Protein
Catalog #: PHH3882

