


Product Name: Recombinant Human KITLG Protein
Catalog #: PHH3908



Summary

Name	SCF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21583
Host	HEK293
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDLKKSFKSP ERLFTPEEF FRIFNRSIDA FKDFVVASET SDCVVSSTLS PEKDSRVSVT KPFMLPPVAA SSLRNDSSSS NRKAKNPPGD

image



Background

Alternative Names Kit ligand (Mast cell growth factor) (MGF) (Stem cell factor) (SCF) (c-Kit ligand) [Cleaved into: Soluble KIT ligand (sKITLG)]

Note

For Research Use Only , Not for Diagnostic Use