

Product Name: Recombinant Human FGA Protein
Catalog #: PHH3911





Summary

Name	FIBA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P02671
Host	HEK293
Species	Human
Molecular Mass	133 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPRVVERHQS ACKDSDWPFC SDEDWNYKCP SGRMKGLID EVNQDFTNRI NKLKNSLFEY QKNNKDSHSL TTNIMEILRG DFSSANNRDN TYNRVSEDLR SRIEVLKRRV IEKVQHIQLL QKNVRAQLVD MKRLEVDIDI KIRSCRGSCS RALAREVDLK DYEDQQKQLE QVIAKDLLPS RDRQHLPLIK MKPVPDLVPG NFKSQLQKVP PEWKALTDMP QMRMELERPG GNEITRGGST SYGTGSETES PRNPSSAGSW NSGSSGPGST GNRNPGSSGT GGTATWKPGS SGPGSTGSWN SGSSGTGSTG NQNP GS PRPG STGTWNP GS ERGSAGHWTS ESSVSGSTGQ WHSESGSFRP DSPGSGNARP NNPDWGT FEE VSGNVSPGTR REYHTEKLV SKGDKELRTG KEKVTSGSTT TTRRSCSKTV TKT VIGPDGH KEVTKEVTS EDGSDCPEAM DLGTL SGIGT LDGFRHRHPD EAAFFDTAST GKTFPGFFSP MLGEFVSETE SRGSESGIFT NTKESSHHP GIAEFPSRGK SSSYSKQFTS STSYNRGDST FESKSYKMAD EAGSEADHEG THSTKRGHAK SRPVRDCDDV LQTHPSGTQS GIFNIKLP GS SKIFSVYCDQ ETSLGGWLLI QQRMDGSLNF NRTWQDYKRG FGSLNDEGEG EFWLGNDYLH LLTQRGSVLR VELEDWAGNE AYA EYHFRVG SEAEGYALQV SSYEGTAGDA LIEGSVEEGA EYTSHNNMQF STFDRDADQW EENCAEVYGG GWWYN NCQAA NLNGIYYPGG SYDPRNNSPY EIENG VVWVS FRGADYSLRA VRMKIRPLVT Q

image

Product Name: Recombinant Human FGA Protein
Catalog #: PHH3911



Background

Alternative Names

Fibrinogen alpha chain [Cleaved into: Fibrinopeptide A; Fibrinogen alpha chain]

Note

For Research Use Only , Not for Diagnostic Use