



Summary

Name	FXYD5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q96DB9
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQTQTQQLEG TDGPLVTDPE THKSTKAAHP TDDTTTLSEK PPSSTDVQTD PQTLKPSGFH EDDPFFYDEH

image



Background

Alternative Names FXYD domain-containing ion transport regulator 5 (Dysadherin)

Note

For Research Use Only , Not for Diagnostic Use