



Product Name: Recombinant Human TNFRSF10B Protein
Catalog #: PHH3944



Summary

Name	TR10B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14763
Host	HEK293
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGEVELSPCT TTRNTVCQCE EGTFREEDSP EMCRCRTGC PRGMVKVGDC TPWSDIECVH KESGTKHSGE

image



Background

Alternative Names	Tumor necrosis factor receptor superfamily member 10B (Death receptor 5) (TNF-related apoptosis-inducing ligand receptor 2) (TRAIL receptor 2) (TRAIL-R2) (CD antigen CD262)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use