



Product Name: Recombinant Human LEPR Protein
Catalog #: PHH3947



Summary

Name	LEPR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P48357
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGHYETAVEP KFNSSGTHFS NLSKTTFHCC FRSEQDRNCS LCADNIEGKT FVSTVNSLVF QQIDANWNIQ CWLKGDLKLF ICYVESLFKN LFRNYNYKVH LLYVLPEVLE DSPLVPQKGS

image



Background

Alternative Names Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295)

Note

For Research Use Only , Not for Diagnostic Use