

Product Name: Recombinant Human L1CAM Protein
Catalog #: PHH3974



Summary

Name	L1CAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P32004
Host	HEK293
Species	Human
Molecular Mass	86 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPRRLVFPT DDISLKCEAS GKPEVQFRWT RDGVHFKPKE ELGVTVYQSP HSGSFTITGN NSNFAQRFGG IYRCFASNKL GTAMSHEIRL MAEGAPKWPK ETVKPVEVEE GESVVLPCNP PPSAEPLRIY WMNSKILHIK QDERVTMGQN GNLYFANVLT SDNHSDYICH AHFPGTRTII QKEPIDLRVK ATNSMIDRKP RLLFPTNSSS HLVALQGQPL VLECIAEGFP TPTIKWLRPS GPMPPADRVTY QNHNKTLQLL KVGEEDDGEY RCLAENSLGS ARHAYYVTV E AAPYWLHKPQ SHLYGPGETA RLDCQVQGRP QPEVTWRING IPVEELAKDQ KYRIQRGALI LSNVQPSDTM VTQCEARNRH GLLANAYIY

image



Background

Alternative Names	Neural cell adhesion molecule L1 (N-CAM-L1) (NCAM-L1) (CD antigen CD171)
--------------------------	--

Product Name: Recombinant Human L1CAM Protein
Catalog #: PHH3974



Note

For Research Use Only , Not for Diagnostic Use