

Product Name: Recombinant Human L1CAM Protein
Catalog #: PHH3975





Summary

Name	L1CAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32004
Host	HEK293
Species	Human
Molecular Mass	124 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLCSPCLLI QIPEEYEGHH VMEPPVITEQ SPRRLVVFPT DDISLKCEAS GKPEVQFRWT RDGVHFKPKE ELGVTVYQSP HSGSFTITGN NSNFAQRFQG IYRCFASNKL GTAMSHEIRL MAEGAPKWPK ETVKPVEVEE GESVVLPCNP PPSAEPLRIY WMNSKILHIK QDERVTMGQN GNLYFANVLT SDNHSDYICH AHFPGTRTII QKEPIDLRVK ATNSMIDRKP RLLFPTNSSS HLVALQGQPL VLECIAEGFP TPTIKWLRPS GPMPADRVTY QNHNKTLLQLL KVGEEDDGEY RCLAENSLGS ARHAYYVTVE AAPYWLHKPQ SHLYGPGETA RLDCQVQGRP QPEVTWRING IPVEELAKDQ KYRIQRGALI LSNVQPSDTM VTQCEARNRH GLLLANAYIY VVQLPAKILT ADNQTYMAVQ GSTAYLLCKA FGAPVPSVQW LDEDGTTVLQ DERFFPYANG TLGIRDQLAN DTGRYFCLAA NDQNNVTIMA NLKVKDATQI TQGPRSTIEK KGSRTVFTCQ ASFDPSLQPS ITWRGDGRDL QELGSDSKYF IEDGRLVIHS LDYSDQGNYS CVA TELDVV ESRAQLLVVG SPGPVPRVLV SDLHLLTQSQ VRVSWSPAED HNAPIEKYDI EFEDKEMAPE KWYSLGKVPV NQTSTTLKLS PYVHYTFRVT AINKYGPGEF SPVSETVVTP EAAPEKNPVD VKGEGNETTN MVITWKPLRW MDWNAPQVQY RVQWRPQGR GPWQEQIVSD PFLVVSNTST FVPYEIKVQA VNSQGKGPEP QVTIGYSGED YPQAIPELEG IEILNSSAVL VKWRPVDLAQ VKGHLRGYNV TYWREGSQRK HSKRHIHKDH VVVPANTTSV ILSGLRPYSS YHLEVQAFNG RGSGPASEFT FSTPEGVPGH PEALHLECQS NTSLLLRWQP PLSHNGVLTG YVLSYHPLDE GGKGQLSFNL RDPELRTHNL TDLSPHLRYR FQLQATTKEG PGEAIVREGG TMALSGISDF GNISATAGEN YSVVSWVPKE GQCNRFRHIL FKALGEEKGG ASLSPQYVSY NQSSYTQWDL QPDTDYEIHL FKERMFRHQM AVKTNGTGRV RLPPAGFATE

image

Product Name: Recombinant Human L1CAM Protein
Catalog #: PHH3975



Background

Alternative Names

Neural cell adhesion molecule L1 (N-CAM-L1) (NCAM-L1) (CD antigen CD171)

Note

For Research Use Only , Not for Diagnostic Use