

Product Name: Recombinant Human CDH16 Protein
Catalog #: PHH4006




Summary

Name	CAD16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O75309
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTVSIPELSP PGTEVTRLA EDADAPGSPN SHVYQLLSP EPEDGVEGRA FQVDPTSGSV TLGVLPLRAG QNILLVLAM DLAGAEGGFS STCEVEVAVT DINDHAPEFI

image



SDS



WB

Background

Alternative Names Cadherin-16 (Kidney-specific cadherin) (Ksp-cadherin)

Note

For Research Use Only , Not for Diagnostic Use