


Product Name: Recombinant Human BMPR1A Protein
Catalog #: PHH4046



Summary

| | |
|------------------------------|---|
| Name | BMR1A_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P36894 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 10 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | QKKSENGVTL APEDTLPLK CYCSGHCPDD AINNTCITNG HCFAIIEEDD QGETTLASGC MKYEGSDFQC |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Bone morphogenetic protein receptor type-1A (BMP type-1A receptor) (BMPR-1A) (EC 2.7.11.30) (Activin receptor-like kinase 3) (ALK-3) (Serine/threonine-protein kinase receptor R5) (SKR5) (CD antigen CD292) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use