


Product Name: Recombinant Human LRP11 Protein
Catalog #: PHH4067



Summary

Name	LRP11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86VZ4
Host	HEK293
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAAPLSELHA QLSGVEQLLE EFRRQLQQR PQEELELELR AGGGPQEDCP GPGSGGYSAM PDAIIRTKDS LAAGASFLRA PAAVRGWRQC VAACCSEPRC SVAVVELPRR PAPPAAVLGC YLFNCTARGR NVCKFALHSG YSSYSLSRAP DGAALATARA SPRQEKDAPP LSKAGQDVVL HLPTDGVVLD GRESTDHAI VQYEWALLQG DPSVDMKVPQ SGTCLKSLHQ EGTYTFQLTV TDTAGQRSSD NVSVTVLRAA YSTGGCLHTC SRYHFFCDDG CCIDITLACD GVQQCPDGSD EDFCQNLGLD RKMVTHTAAS PALPRTTGPS EDAGGDSLVE KSQKATAPNK PPALSNTKR NHSAFWGPES QIIPVMPDSS SSGKNRKEES YIFESKGDGG GGEHPAPETG

image



Background

Product Name: Recombinant Human LRP11 Protein
Catalog #: PHH4067



Alternative Names

Low-density lipoprotein receptor-related protein 11 (LRP-11)

Note

For Research Use Only , Not for Diagnostic Use