


Product Name: Recombinant Human CXCL16 Protein
Catalog #: PHH4084



Summary

Name	CXL16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2A7
Host	HEK293
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGSVTGSCYC GKRISDSPP SVQFMNRLRK HLRAYHRCLY YTRFQLLSWS VCGGNKDPWV QELMSCLDLK ECGHAYSGIV AHQKHLLPTS PPISQASEGA SSDIHTPAQM LLSTLQSTQR PTLPVGSLSS DKELTRPNET TIHTAGHSLA AGPEAGENQK QPEKNAGPTA

image



Background

Alternative Names	C-X-C motif chemokine 16 (Scavenger receptor for phosphatidylserine and oxidized low density lipoprotein) (SR-PSOX) (Small-inducible cytokine B16) (Transmembrane chemokine CXCL16)
--------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

Product Name: Recombinant Human CXCL16 Protein
Catalog #: PHH4084



For Research Use Only , Not for Diagnostic Use